



# nSequence<sup>®</sup> LITE

## A Streamlined Option for Your Full Arch Restorations

**When immediate loading may not be an option, nSequence Lite brings an easier, faster and more predictable solution.**

- » Enhances surgical accuracy while ensuring proper bone reduction for a Definitive Prosthesis.
- » Reduces restorative appointments for the Definitive Prosthesis and provides a streamlined workflow.
- » The Clear Duplicate Denture functions as an all-in-one verification jig, bite block, set-up/try-in and final impression.
- » Opens the door for many overdenture options, including Locators, Hader Bar<sup>®</sup>, Conus, and for hybrid Definitive Restorations.

**ALL YOU NEED ARE 4 COMPONENTS:**

1 BONE FOUNDATION GUIDE WITH MONOSTRUT	2 SURGICAL GUIDE	3 CLEAR DUPLICATE DENTURE	4 IMMEDIATE DENTURE
			

**SEND US YOUR CASE TODAY! [www.nSequence.com](http://www.nSequence.com)**